

Table of Contents

General Information

Notice to our Customers	2-3
-----------------------------------	-----

Prosthetic Products

General Product Information	4
How to Size Royal Knit Prosthetic Socks	5
How to Size Royal Knit Symes Socks	5
Royal Knit's Easy to Use Wall Chart	6
Super Sheath (Heavy and Lightweight)	7-8
Underliner Sock	9
Double Fleece Socks	10
The Royal Flush Sock - 3, 5, 6 & 8 ply	11-12
The Royal 12 Pack	13
The Royal Open Toe	14
The Royal Flush Multi-Ply 3, 5, 6 Ply Combination Sock	15
The Royal Stretch Sock - 3, 5 and 6 ply	16-17
Post OP Sock	18
Hemipelvectomy/Hip Disarticulation Sock	19

Orthotic Products

Torso Body Sock	20-21
Fracture Interface Brace/Cast Socks	22
Full Fashioned Fracture/Brace Sock	23
Diabetic Sock	23
Cranial Caps (1 & 2 Ply)	23
Fracture Interface Arm Sleeve, Humeral & Ulnar Sleeves	24
Hip Orthosis Sock	25
Toe Sock	25

L-Codes	26
--------------------------	----

Notes	27
------------------------	----

NOTICE TO OUR CUSTOMERS

QUICK RESPONSE

Royal Knit is known throughout our industry for the standard of excellence we have established in immediately filling customer orders. Items in stock will be shipped the same day if your order is in before 3:30 p.m. (CST). Royal Knit has ninety-nine percent (99%) of their items available for shipment on any given day.

ORDERING PROCEDURE

There are three ways to order Royal Knit socks. Orders may be placed by telephone, fax or e-mail royalknitinc@yahoo.com. To expedite order processing for socks please indicate the following:

Your account number, material, ply, size group, length and quantity.

ORDERING EXAMPLES

RF3-0/1-10 (Quantity 12)

RF3 = Royal Flush 3 ply

0/1 = Size (small, 6"-7" top, 3"-4" toe)

10 = Length of 8-10"

Quantity 12

RS3-0/1-10 (Quantity 12)

RS3 = Stretch Wool 3 ply

EXAMPLE WITH DISTAL OPENING (SAME SOCK)

Please notice ordering a sock with a distal opening, simply requires adding an "O" before the product code.

EXAMPLE WITH RUBBERIZED REINFORCED DISTAL OPENING WITH (SAME SOCK)

If you want a rubberized, reinforced distal opening add a "GO" before the product code.

TERMS

Standard credit terms are thirty (30) days unless specifically stated otherwise on invoice. For accounts which are past due, Royal Knit reserves the right to assess a maximum charge of one and one-half percent (1.5% per month) on the outstanding balance. (18% annually)

MASTERCARD/VISA PURCHASES

Customers may place orders using Mastercard or Visa. Please notify your Royal Knit representative when ordering by phone. **Credit card orders do not qualify for early pay discounts.**

PRICES

All prices are in U.S. dollars. Prices are subject to change without notice. All orders are FOB shipping point. If an item is in stock it will be invoiced the next business day. If an item is not in stock it will be invoiced upon shipment.

SPECIAL ORDERS

Items ordered as special orders or items not in stock are non-cancelable and non-returnable. Special orders and items not in stock are normally shipped in one (1) business days.

DOMESTIC SHIPPING

All orders are automatically shipped UPS Ground unless otherwise specified.

RETURNS

There is a strict 30 day return policy. All returns must have a return authorization (RA) number and we have to receive the items back within thirty-days (30) of the invoice date. ***THE PRODUCT MUST BE IN IT'S ORIGINAL PACKAGE AND MUST NOT HAVE BEEN OPENED.*** Contact your Royal Knit representative for an RA number. Further instructions may be given at that time.

NON-RETURNABLE ITEMS INCLUDE: (ALL LISTED ARE MADE TO ORDER)

- 8 ply socks
- Split ply socks
- Socks 22" and above in length
- Custom special order
- Royal Silver
- Modified products (open toe socks for seal-in liners, sewn tops, socks cut to length)

If you aren't sure, please ask your customer service representative if your products are returnable.

GENERAL PRODUCT INFORMATION

The product information included below is designed to explain how our prosthetic socks are sized, the advantages in our sizing method, and the fiscal benefits in carrying the Royal Knit line.

STRETCH

When sizing our socks, taking the degree of stretch into consideration is paramount. A stretch sock, much like a rubber band, gets thinner when stretched. This may cause a loss of ply. The amount of stretch will determine the amount of ply lost and will vary from patient to patient. Understanding this, Royal Knit offers our prosthetic sock lines in overlapping sizing.

OVERLAPPING SIZING

How Overlapping Sizing Works

Due to our overlapping sizes, we initially recommend trying the smallest size possible when fitting the patient. Knowing that stretch allows for a more conforming fit, it is natural that a smaller size that fits will provide a snugger, more comfortable sock.

For example, if your patient normally wears #2-12, we suggest our size 1/2-12 which allows a three inch stretch. If the patient feels there is substantial loss of ply then we suggest trying size 2/3-12. This will allow a one and one-half inch stretch circumference. Since less stretch is required in the larger size the loss of ply will be minimized.

The concept of overlapping fit can be taken one step further. If the sock used in the example above is a five ply product and the patient feels the 2/3-12 provides a better fit after trying the 1/2-12, but still feels there is a lack of ply in the sock, it is at this point we suggest moving to a heavy 5- ply. Again providing another opportunity to minimize the loss of ply.

Product Qualities

In addition to the overlapping size concept, there are other advantages in carrying the Royal Knit line.

- Easier application for the geriatric patient.
- Stretch and full fashioned knitting techniques provide a snug, comfortable fit for almost any size. Usually no special sizing required.
- Most standard sizes fit both PTB and SYMES patients.
- Rounded toe (no dog ears).
- Ply stamp that will last the life of the product.
- Reinforced Rubberized Hole will last the life of the product.
- The stretch in our socks offers the advantage of carrying fewer sizes, thus reducing the inventory necessary to accommodate patient needs.
- Wash and Wear
- Hemipelvectomy/Hip disarticulation Sock are stock items, waist measurements are needed.

SUGGESTED INVENTORY

Stocking Royal Knit prosthetic socks offers unique opportunities in budgetary management.

Due to the natural stretch of our socks we have been able to condense 132 non stretch sizes into 48 overlapping sizes. For example, facilities can reduce inventory by stocking the following:

Recommended core inventory socks.

60% - 1/2-12 BK

25% - 2/3-14 BK/AK

15% - 3/4-14 AK

INVENTORY CONVERSION

If you have any problems in converting our products codes sizing or other information from different sizing methods, Royal Knit will be happy to assist you. Just provide your current inventory. If you elect, Royal Knit will also make recommendations for inventory requirements and potential reductions.

CARE INSTRUCTIONS


There are no special laundering instructions. Simply machine wash with warm water, and tumble dry, using the setting tumble dry warm. This will bring the sock back to the original stated measurements. Socks including Body Socks must be machined washed and dried.

Royal Knit socks (white) **can be bleached**, in cold water. Adding recommended amounts of bleach during washing will not harm the sock.


Royal Knit socks are produced with resistance to shrinking, odors and mildew. If properly cared for, the Royal Flush and wool line of socks will retain their shape, elasticity and wicking action after repeated washing.

Sizing Chart

Standard Sock



Symes Sock



THIS CHART APPLIES ONLY TO 3,5,6,8 PLY STANDARD
ALSO APPLIES TO 3,5,6 HEAVY, AND WOOL

SIZE	TOP	TOE	LENGTHS
A/B = (x-small)	4"-4.5"	3"-3.5"	6", 10", 12", 14", 18", 22"
B/0 = (x-small/sm.)	5"-6"	3"-3.5"	10", 12", 14", 18", 22"
0/1 = (small)	6"-7"	3"-4"	10", 12", 14", 18", 22", 26"
0/2 = (sm./md. ptb)	6"-7"	4"-5"	10", 12", 14", 18", 22", 26"
1/2 = (medium)	7"-8"	4"-5"	10", 12", 14", 18", 22", 26"
2/3 = (large)	8"-9"	5"-6"	10", 12", 14", 18", 22", 26"
3/4 = (x-large)	9"-10"	6"-7"	10", 12", 14", 18", 22", 26"
4/6 = (xx-large)	10"-12"	7"-8"	10", 12", 14", 18", 22", 26"

THIS CHART APPLIES TO SUPER SHEATHS, AND HEAVY SHEATHS

SIZE	TOP	TOE	LENGTHS
A/B = (x-small)	4"-4.5"	3"-3.5"	6", 10", 12", 14", 16", 20"
B/0 = (x-small/sm)	5"-6"	3"-3.5"	10", 12", 14", 16", 20"
0/1 = (small)	6"-7"	3"-4"	10", 12", 14", 16", 20", 24"
0/2 = (sm./md. ptb)	6"-7"	4"-5"	10", 12", 14", 16", 20", 24"
1/2 = (medium)	7"-8"	4"-5"	10", 12", 14", 16", 20", 24"
2/3 = (large)	8"-9"	5"-6"	10", 12", 14", 16", 20", 24"
3/4 = (x-large)	9"-10"	6"-7"	10", 12", 14", 16", 20", 24"
4/6 = xx-large)	10"-12"	7"-8"	10", 12", 14", 16", 20", 24"

THIS CHART APPLIES TO THE 1 AND 2 PLY

SIZE	TOP	TOE	LENGTHS
0A = (x-small)	4"-5"	2.5"-3"	8", 14"
01 = (small)	5.5"-7"	3"-4"	10", 12", 16", 20", 24"
12 = (medium)	7"-8"	4"-5"	10", 12", 16", 20", 28"
23 = (large)	8"-9"	5"-6"	10", 12", 16", 20"
35 = (x-large)	9"-11"	5.5"-7"	10", 12", 16", 20"

ROYAL KNIT'S WALL CHART (3, 5, 6 and 8 ply only)

X-Small				
Product Code	Size	Length	Top	Toe
__RF__-A/B-6	#A	4"-6"	4"-4½"	3"-3½"
__RF__-A/B-10	#A	8"-10"	4"-4½"	3"-3½"
__RF__-A/B-12	#A	10"-12"	4"-4½"	3"-3½"
__RF__-A/B-14	#A	12"-14"	4"-4½"	3"-3½"
__RF__-A/B-18	#A	16"-18"	4"-4½"	3"-3½"
__RF__-A/B-22	#A	20"-22"	4"-4½"	3"-3½"

CHILD/ARM

**LARGE B/K
SMALL A/K**
(fits #2 or size #3)

X-Small/Small				
Product Code	Size	Length	Top	Toe
__RF__-B/0-10	#B-#0	8"-10"	5"-6"	3"-3¾"
__RF__-B/0-12	#B-#0	10"-12"	5"-6"	3"-3¾"
__RF__-B/0-14	#B-#0	12"-14"	5"-6"	3"-3¾"
__RF__-B/0-18	#B-#0	16"-18"	5"-6"	3"-3¾"
__RF__-B/0-22	#B-#0	20"-22"	5"-6"	3"-3¾"

CHILD/ARM

(same toe / larger top)

STANDARD A/K
(fits #3 or size #4)

Small				
Product Code	Size	Length	Top	Toe
__RF__-0/1-10	#0-#1	8"-10"	6"-7"	3"-4"
__RF__-0/1-12	#0-#1	10"-12"	6"-7"	3"-4"
__RF__-0/1-14	#0-#1	12"-14"	6"-7"	3"-4"
__RF__-0/1-18	#0-#1	16"-18"	6"-7"	3"-4"
__RF__-0/1-22	#0-#1	20"-22"	6"-7"	3"-4"
__RF__-0/1-26	#0-#1	24"-26"	6"-7"	3"-4"

SMALL B/K

(fits size #0 or size #1)

EXTREMELY LARGE A/K

Small/Medium (PTB)				
Product Code	Size	Length	Top	Toe
__RF__-0/2-10	Slim	8"-10"	6"-7"	4"-5"
__RF__-0/2-12	Slim	10"-12"	6"-7"	4"-5"
__RF__-0/2-14	Slim	12"-14"	6"-7"	4"-5"
__RF__-0/2-18	Slim	16"-18"	6"-7"	4"-5"
__RF__-0/2-22	Slim	20"-22"	6"-7"	4"-5"
__RF__-0/2-26	Slim	24"-26"	6"-7"	4"-5"

SYMMETRICAL

(size medium toe)
(size small top)



STANDARD B/K

(fits size #1 or size #2)

Medium				
Product Code	Size	Length	Top	Toe
__RF__-1/2-10	#1-#2	8"-10"	7"-8"	4"-5"
__RF__-1/2-12	#1-#2	10"-12"	7"-8"	4"-5"
__RF__-1/2-14	#1-#2	12"-14"	7"-8"	4"-5"
__RF__-1/2-18	#1-#2	16"-18"	7"-8"	4"-5"
__RF__-1/2-22	#1-#2	20"-22"	7"-8"	4"-5"
__RF__-1/2-26	#1-#2	24"-26"	7"-8"	4"-5"

Latex Free

98% Tetra/channel polyester 2% Lycra®

Large				
Product Code	Size	Length	Top	Toe
__RF__-2/3-10	#2-#3	8"-10"	8"-9"	5"-6"
__RF__-2/3-12	#2-#3	10"-12"	8"-9"	5"-6"
__RF__-2/3-14	#2-#3	12"-14"	8"-9"	5"-6"
__RF__-2/3-18	#2-#3	16"-18"	8"-9"	5"-6"
__RF__-2/3-22	#2-#3	20"-22"	8"-9"	5"-6"
__RF__-2/3-26	#2-#3	24"-26"	8"-9"	5"-6"

X-Large				
Product Code	Size	Length	Top	Toe
__RF__-3/4-10	#3-#4	8"-10"	9"-10"	6"-7"
__RF__-3/4-12	#3-#4	10"-12"	9"-10"	6"-7"
__RF__-3/4-14	#3-#4	12"-14"	9"-10"	6"-7"
__RF__-3/4-18	#3-#4	16"-18"	9"-10"	6"-7"
__RF__-3/4-22	#3-#4	20"-22"	9"-10"	6"-7"
__RF__-3/4-26	#3-#4	24"-26"	9"-10"	6"-7"

XX-Large				
Product Code	Size	Length	Top	Toe
__RF__-4/6-10	#5-#6	8"-10"	10"-12"	7"-8"
__RF__-4/6-12	#5-#6	10"-12"	10"-12"	7"-8"
__RF__-4/6-14	#5-#6	12"-14"	10"-12"	7"-8"
__RF__-4/6-18	#5-#6	16"-18"	10"-12"	7"-8"
__RF__-4/6-22	#5-#6	20"-22"	10"-12"	7"-8"
__RF__-4/6-26	#5-#6	24"-26"	10"-12"	7"-8"

CARE INSTRUCTIONS

- Machine wash warm
- Bleach if needed, cold water only
- Pre-treat if needed
- Tumble dry warm
- Use a sheet fabric softener (no liquid softener)
- Do not hand wash

Super Sheath

The Super Sheath works great in elevated vacuum, underliner socks or a thinner fitting sock for liners. It is designed to offer a form fitting prosthetic sheath. This creates a fuller, smoother sheath for increased cushioning and softness. The round toe is engineered for a “knit-to-fit” feeling with consistent memory and excellent stretch. If properly cared for, stays soft and resists shrinking.

L-CODE: L8410 - AK L-CODE: L8400 - BK

See **OPTIONS** below for our easy add-ons

REINFORCED DISTAL OPENING-(Add “GO” in front of product code) Ex: (**GORN**S)
OPEN TOE-(Add “S” in front of product code) Ex: (**SR**NS)

- Stretches in Length and Width
- Form Fitting
- Most Sizes Fit Symes and PTB Patients
- Round Toe
- Wash and Dry Low Heat



A/B (X-SMALL)		4”-5” TOP X 3”-3.5” TOE- FLAT MEASUREMENT							
NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE	
RNS-A/B-06	ORNS-A/B-06	BRNS-A/B-06	OBRNS-A/B-06	HRNS-A/B-06	OHRNS-A/B-06	BHRNS-A/B-06	OBHRNS-A/B-06	CHILD X-SHORT	
RNS-A/B-10	ORNS-A/B-10	BRNS-A/B-10	OBRNS-A/B-10	HRNS-A/B-10	OHRNS-A/B-10	BHRNS-A/B-10	OBHRNS-A/B-10	CHILD SHORT	
RNS-A/B-12	ORNS-A/B-12	BRNS-A/B-12	OBRNS-A/B-12	HRNS-A/B-12	OHRNS-A/B-12	BHRNS-A/B-12	OBHRNS-A/B-12	CHILD MEDIUM	
RNS-A/B-14	ORNS-A/B-14	BRNS-A/B-14	OBRNS-A/B-14	HRNS-A/B-14	OHRNS-A/B-14	BHRNS-A/B-14	OBHRNS-A/B-14	RK ORIGINAL	
RNS-A/B-16	ORNS-A/B-16	BRNS-A/B-16	OBRNS-A/B-16	HRNS-A/B-16	OHRNS-A/B-16	BHRNS-A/B-16	OBHRNS-A/B-16	RK ORIGINAL	
RNS-A/B-20	ORNS-A/B-20	BRNS-A/B-20	OBRNS-A/B-20	HRNS-A/B-20	OHRNS-A/B-20	BHRNS-A/B-20	OBHRNS-A/B-20	RK ORIGINAL	

B/O (X-SMALL / SMALL)		5”-6” TOP X 3”-3.5” TOE- FLAT MEASUREMENT							
NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE	
RNS-B/0-10	ORNS-B/0-10	BRNS-B/0-10	OBRNS-B/0-10	HRNS-B/0-10	OHRNS-B/0-10	BHRNS-B/0-10	OBHRNS-B/0-10	RK ORIGINAL	
RNS-B/0-12	ORNS-B/0-12	BRNS-B/0-12	OBRNS-B/0-12	HRNS-B/0-12	OHRNS-B/0-12	BHRNS-B/0-12	OBHRNS-B/0-12	RK ORIGINAL	
RNS-B/0-14	ORNS-B/0-14	BRNS-B/0-14	OBRNS-B/0-14	HRNS-B/0-14	OHRNS-B/0-14	BHRNS-B/0-14	OBHRNS-B/0-14	RK ORIGINAL	
RNS-B/0-16	ORNS-B/0-16	BRNS-B/0-16	OBRNS-B/0-16	HRNS-B/0-16	OHRNS-B/0-16	BHRNS-B/0-16	OBHRNS-B/0-16	RK ORIGINAL	
RNS-B/0-20	ORNS-B/0-20	BRNS-B/0-20	OBRNS-B/0-20	HRNS-B/0-20	OHRNS-B/0-20	BHRNS-B/0-20	OBHRNS-B/0-20	RK ORIGINAL	

O/1 (SMALL)		6”-7” TOP X 3”-4” TOE- FLAT MEASUREMENT							
NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE	
RNS-0/1-10	ORNS-0/1-10	BRNS-0/1-10	OBRNS-0/1-10	HRNS-0/1-10	OHRNS-0/1-10	BHRNS-0/1-10	OBHRNS-0/1-10	NARROW X-SHORT	
RNS-0/1-12	ORNS-0/1-12	BRNS-0/1-12	OBRNS-0/1-12	HRNS-0/1-12	OHRNS-0/1-12	BHRNS-0/1-12	OBHRNS-0/1-12	NARROW SHORT	
RNS-0/1-14	ORNS-0/1-14	BRNS-0/1-14	OBRNS-0/1-14	HRNS-0/1-14	OHRNS-0/1-14	BHRNS-0/1-14	OBHRNS-0/1-14	NARROW SHORT	
RNS-0/1-16	ORNS-0/1-16	BRNS-0/1-16	OBRNS-0/1-16	HRNS-0/1-16	OHRNS-0/1-16	BHRNS-0/1-16	OBHRNS-0/1-16	NARROW MEDIUM	
RNS-0/1-20	ORNS-0/1-20	BRNS-0/1-20	OBRNS-0/1-20	HRNS-0/1-20	OHRNS-0/1-20	BHRNS-0/1-20	OBHRNS-0/1-20	NARROW MEDIUM	
RNS-0/1-24	ORNS-0/1-24	BRNS-0/1-24	OBRNS-0/1-24	HRNS-0/1-24	OHRNS-0/1-24	BHRNS-0/1-24	OBHRNS-0/1-24	NARROW LONG	

0/2 (SMALL/MEDIUM)**6"-7" TOP X 4"-5" TOE- FLAT MEASUREMENT**

NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE
RNS-0/2-10	ORNS-0/2-10	BRNS-0/2-10	OBRNS-0/2-10	HRNS-0/2-10	OHRNS-0/2-10	BHRNS-0/2-10	OBHRNS-0/2-10	RK ORIGINAL
RNS-0/2-12	ORNS-0/2-12	BRNS-0/2-12	OBRNS-0/2-12	HRNS-0/2-12	OHRNS-0/2-12	BHRNS-0/2-12	OBHRNS-0/2-12	RK ORIGINAL
RNS-0/2-14	ORNS-0/2-14	BRNS-0/2-14	OBRNS-0/2-14	HRNS-0/2-14	OHRNS-0/2-14	BHRNS-0/2-14	OBHRNS-0/2-14	RK ORIGINAL
RNS-0/2-16	ORNS-0/2-16	BRNS-0/2-16	OBRNS-0/2-16	HRNS-0/2-16	OHRNS-0/2-16	BHRNS-0/2-16	OBHRNS-0/2-16	RK ORIGINAL
RNS-0/2-20	ORNS-0/2-20	BRNS-0/2-20	OBRNS-0/2-20	HRNS-0/2-20	OHRNS-0/2-20	BHRNS-0/2-20	OBHRNS-0/2-20	RK ORIGINAL
RNS-0/2-24	ORNS-0/2-24	BRNS-0/2-24	OBRNS-0/2-24	HRNS-0/2-24	OHRNS-0/2-24	BHRNS-0/2-24	OBHRNS-0/2-24	RK ORIGINAL

1/2 (MEDIUM)**7"-8" TOP X 4"-5" TOE- FLAT MEASUREMENT**

NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE
RNS-1/2-10	ORNS-1/2-10	BRNS-1/2-10	OBRNS-1/2-10	HRNS-1/2-10	OHRNS-1/2-10	BHRNS-1/2-10	OBHRNS-1/2-10	REGULAR X-SHORT
RNS-1/2-12	ORNS-1/2-12	BRNS-1/2-12	OBRNS-1/2-12	HRNS-1/2-12	OHRNS-1/2-12	BHRNS-1/2-12	OBHRNS-1/2-12	REGULAR SHORT
RNS-1/2-14	ORNS-1/2-14	BRNS-1/2-14	OBRNS-1/2-14	HRNS-1/2-14	OHRNS-1/2-14	BHRNS-1/2-14	OBHRNS-1/2-14	REGULAR SHORT
RNS-1/2-16	ORNS-1/2-16	BRNS-1/2-16	OBRNS-1/2-16	HRNS-1/2-16	OHRNS-1/2-16	BHRNS-1/2-16	OBHRNS-1/2-16	REGULAR MEDIUM
RNS-1/2-20	ORNS-1/2-20	BRNS-1/2-20	OBRNS-1/2-20	HRNS-1/2-20	OHRNS-1/2-20	BHRNS-1/2-20	OBHRNS-1/2-20	REGULAR MEDIUM
RNS-1/2-24	ORNS-1/2-24	BRNS-1/2-24	OBRNS-1/2-24	HRNS-1/2-24	OHRNS-1/2-24	BHRNS-1/2-24	OBHRNS-1/2-24	REGULAR LONG

2/3 (LARGE)**8"-9" TOP X 5"-6" TOE- FLAT MEASUREMENT**

NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE
RNS-2/3-10	ORNS-2/3-10	BRNS-2/3-10	OBRNS-2/3-10	HRNS-2/3-10	OHRNS-2/3-10	BHRNS-2/3-10	OBHRNS-2/3-10	RK ORIGINAL
RNS-2/3-12	ORNS-2/3-12	BRNS-2/3-12	OBRNS-2/3-12	HRNS-2/3-12	OHRNS-2/3-12	BHRNS-2/3-12	OBHRNS-2/3-12	RK ORIGINAL
RNS-2/3-14	ORNS-2/3-14	BRNS-2/3-14	OBRNS-2/3-14	HRNS-2/3-14	OHRNS-2/3-14	BHRNS-2/3-14	OBHRNS-2/3-14	RK ORIGINAL
RNS-2/3-16	ORNS-2/3-16	BRNS-2/3-16	OBRNS-2/3-16	HRNS-2/3-16	OHRNS-2/3-16	BHRNS-2/3-16	OBHRNS-2/3-16	RK ORIGINAL
RNS-2/3-20	ORNS-2/3-20	BRNS-2/3-20	OBRNS-2/3-20	HRNS-2/3-20	OHRNS-2/3-20	BHRNS-2/3-20	OBHRNS-2/3-20	RK ORIGINAL
RNS-2/3-24	ORNS-2/3-24	BRNS-2/3-24	OBRNS-2/3-24	HRNS-2/3-24	OHRNS-2/3-24	BHRNS-2/3-24	OBHRNS-2/3-24	RK ORIGINAL

3/4 (X-LARGE)**9"-10" TOP X 6"-7" TOE- FLAT MEASUREMENT**

NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE
RNS-3/4-10	ORNS-3/4-10	BRNS-3/4-10	OBRNS-3/4-10	HRNS-3/4-10	OHRNS-3/4-10	BHRNS-3/4-10	OBHRNS-3/4-10	WIDE X-SHORT
RNS-3/4-12	ORNS-3/4-12	BRNS-3/4-12	OBRNS-3/4-12	HRNS-3/4-12	OHRNS-3/4-12	BHRNS-3/4-12	OBHRNS-3/4-12	WIDE SHORT
RNS-3/4-14	ORNS-3/4-14	BRNS-3/4-14	OBRNS-3/4-14	HRNS-3/4-14	OHRNS-3/4-14	BHRNS-3/4-14	OBHRNS-3/4-14	WIDE SHORT
RNS-3/4-16	ORNS-3/4-16	BRNS-3/4-16	OBRNS-3/4-16	HRNS-3/4-16	OHRNS-3/4-16	BHRNS-3/4-16	OBHRNS-3/4-16	WIDE MEDIUM
RNS-3/4-20	ORNS-3/4-20	BRNS-3/4-20	OBRNS-3/4-20	HRNS-3/4-20	OHRNS-3/4-20	BHRNS-3/4-20	OBHRNS-3/4-20	WIDE MEDIUM
RNS-3/4-24	ORNS-3/4-24	BRNS-3/4-24	OBRNS-3/4-24	HRNS-3/4-24	OHRNS-3/4-24	BHRNS-3/4-24	OBHRNS-3/4-24	RK ORIGINAL

4/6 (XX-LARGE)**10"-12" TOP X 7"-8" TOE- FLAT MEASUREMENT**

NON-HOLE WHITE	HOLE-IN-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	NON-HOLE WHITE/HEAVY	HOLE-IN-TOE WHITE/HEAVY	NON-HOLE BLACK/HEAVY	HOLE-IN-TOE BLACK/HEAVY	COMP SIZE
RNS-4/6-10	ORNS-4/6-10	BRNS-4/6-10	OBRNS-4/6-10	HRNS-4/6-10	OHRNS-4/6-10	BHRNS-4/6-10	OBHRNS-4/6-10	X-WIDE X-SHORT
RNS-4/6-12	ORNS-4/6-12	BRNS-4/6-12	OBRNS-4/6-12	HRNS-4/6-12	OHRNS-4/6-12	BHRNS-4/6-12	OBHRNS-4/6-12	X-WIDE SHORT
RNS-4/6-14	ORNS-4/6-14	BRNS-4/6-14	OBRNS-4/6-14	HRNS-4/6-14	OHRNS-4/6-14	BHRNS-4/6-14	OBHRNS-4/6-14	X-WIDE SHORT
RNS-4/6-16	ORNS-4/6-16	BRNS-4/6-16	OBRNS-4/6-16	HRNS-4/6-16	OHRNS-4/6-16	BHRNS-4/6-16	OBHRNS-4/6-16	X-WIDE MEDIUM
RNS-4/6-20	ORNS-4/6-20	BRNS-4/6-20	OBRNS-4/6-20	HRNS-4/6-20	OHRNS-4/6-20	BHRNS-4/6-20	OBHRNS-4/6-20	X-WIDE MEDIUM
RNS-4/6-24	ORNS-4/6-24	BRNS-4/6-24	OBRNS-4/6-24	HRNS-4/6-24	OHRNS-4/6-24	BHRNS-4/6-24	OBHRNS-4/6-24	RK ORIGINAL

THE ROYAL Liner Sock

The Royal Liner Sock is the perfect under-liner sock for suspension liners. It is to be worn under the liner to relieve any skin irritation and improves the comfort of the liner. The fibers are made to keep the skin dry and cool throughout the day, as well as to eliminate odors.

Features of The Royal Liner Sock:

- **Non-Roll Top!!** - Eliminates bunching for a smoother, more comfortable fit.
- Available in White, Black, Silver.
- Form fitting.
- Wicks moisture away from the skin.
- Easy care and shrinkage resistant.
- Silver fiber helps prevent odors and is antimicrobial.



SIZING CHART

White	Black	Silver	Size	Fit Sock Size	Length
RK1P-01-10	BRK1P-01-10	XRK1P-01-10	Small X-Short	B, 0	8" - 12" (20-31cm)
RK1P-01-12	BRK1P-01-12	XRK1P-01-12	Small Short	B, 0	10" - 14" (25-35cm)
RK1P-12-08	BRK1P-12-08	XRK1P-12-08	Medium XX-Short	1, 2	6" - 10" (15-25cm)
RK1P-12-10	BRK1P-12-10	XRK1P-12-10	Medium X-Short	1, 2	8" - 12" (20-31cm)
RK1P-12-12	BRK1P-12-12	XRK1P-12-12	Medium Short	1, 2	10" - 14" (25-35cm)
RK1P-12-16	BRK1P12-16	XRK1P-12-16	Medium Medium	1, 2	16" - 20" (40-50cm)
RK1P-35-12	BRK1P-35-12	XRK1P-35-12	X-Large Short	3+	10" - 14" (25-35cm)
RK1P-35-16	BRK1P-35-16	XRK1P-35-16	X-Large Medium	3+	16" - 20" (40-50cm)

Other sizes are available.

The Royal Double Fleece

The Royal Double Fleece is designed to quickly move moisture away from the skin to the outer layer of the fabric. The breathability of the fiber keeps you drier and more comfortable. Available in either a 1 and 2 ply.

L-CODE: L8430 - AK L-CODE: L8420 - BK

See **OPTIONS** below for our easy add-ons

SILVER-(Add "X" in front of product code) Ex: (**XRKDF**)

REINFORCED DISTAL OPENING-(Add "GO" in front of product code) Ex: (**GORK1P**)

- Non-Roll Top
- Reinforced Toe with Tighter Stitch Pattern
- Most Sizes Fit Symes and PTB Patients

- Round Toe
- Wash and Dry Low Heat
- Form Fitting



0A (X-SMALL)		4"-5" TOP X 2.5"-3" TOE- FLAT MEASUREMENT							
1 PLY NON-HOLE WHITE	1 PLY HOLE-IN-TOE WHITE	1 PLY NON-HOLE BLACK	1 PLY HOLE-IN-TOE BLACK	2 PLY NON-HOLE WHITE	2 PLY HOLE-IN-TOE WHITE	2 PLY NON-HOLE BLACK	2 PLY HOLE-IN-TOE BLACK	COMP SIZE	
RK1P-0A-08	ORK1P-0A-08	BRK1P-0A-08	OBRK1P-0A-08	RKDF-0A-08	ORKDF-0A-08	BRKDF-0A-08	OBRKDF-0A-08	CHILD SHORT	
RK1P-0A-14	ORK1P-0A-14	BRK1P-0A-14	OBRK1P-0A-14	RKDF-0A-14	ORKDF-0A-14	BRKDF-0A-14	OBRKDF-0A-14	RK ORIGINAL	

01 (SMALL)		5.5"-7" TOP X 3"-4" TOE- FLAT MEASUREMENT							
1 PLY NON-HOLE WHITE	1 PLY HOLE-IN-TOE WHITE	1 PLY NON-HOLE BLACK	1 PLY HOLE-IN-TOE BLACK	2 PLY NON-HOLE WHITE	2 PLY HOLE-IN-TOE WHITE	2 PLY NON-HOLE BLACK	2 PLY HOLE-IN-TOE BLACK	COMP SIZE	
RK1P-01-10	ORK1P-01-10	BRK1P-01-10	OBRK1P-01-10	RKDF-01-10	ORKDF-01-10	BRKDF-01-10	OBRKDF-01-10	NARROW X-SHORT	
RK1P-01-12	ORK1P-01-12	BRK1P-01-12	OBRK1P-01-12	RKDF-01-12	ORKDF-01-12	BRKDF-01-12	OBRKDF-01-12	NARROW SHORT	
RK1P-01-16	ORK1P-01-16	BRK1P-01-16	OBRK1P-01-16	RKDF-01-16	ORKDF-01-16	BRKDF-01-16	OBRKDF-01-16	NARROW MEDIUM	
RK1P-01-20	ORK1P-01-20	BRK1P-01-20	OBRK1P-01-20	RKDF-01-20	ORKDF-01-20	BRKDF-01-20	OBRKDF-01-20	NARROW LONG	
RK1P-01-24	ORK1P-01-24	BRK1P-01-24	OBRK1P-01-24	RKDF-01-24	ORKDF-01-24	BRKDF-01-24	OBRKDF-01-24	NARROW X-LONG	

12 (MEDIUM)		7"-8" TOP X 4"-5" TOE- FLAT MEASUREMENT							
1 PLY NON-HOLE WHITE	1 PLY HOLE-IN-TOE WHITE	1 PLY NON-HOLE BLACK	1 PLY HOLE-IN-TOE BLACK	2 PLY NON-HOLE WHITE	2 PLY HOLE-IN-TOE WHITE	2 PLY NON-HOLE BLACK	2 PLY HOLE-IN-TOE BLACK	COMP SIZE	
RK1P-12-10	ORK1P-12-10	BRK1P-12-10	OBRK1P-12-10	RKDF-12-10	ORKDF-12-10	BRKDF-12-10	OBRKDF-12-10	REGULAR X-SHORT	
RK1P-12-12	ORK1P-12-12	BRK1P-12-12	OBRK1P-12-12	RKDF-12-12	ORKDF-12-12	BRKDF-12-12	OBRKDF-12-12	REGULAR SHORT	
RK1P-12-16	ORK1P-12-16	BRK1P-12-16	OBRK1P-12-16	RKDF-12-16	ORKDF-12-16	BRKDF-12-16	OBRKDF-12-16	REGULAR MEDIUM	
RK1P-12-20	ORK1P-12-20	BRK1P-12-20	OBRK1P-12-20	RKDF-12-20	ORKDF-12-20	BRKDF-12-20	OBRKDF-12-20	REGULAR LONG	
RK1P-12-28	ORK1P-12-28	BRK1P-12-28	OBRK1P-12-28	RKDF-12-28	ORKDF-12-28	BRKDF-12-28	OBRKDF-12-28	REGULAR X-LONG	

23 (LARGE)		8"-9" TOP X 5"-6" TOE- FLAT MEASUREMENT							
1 PLY NON-HOLE WHITE	1 PLY HOLE-IN-TOE WHITE	1 PLY NON-HOLE BLACK	1 PLY HOLE-IN-TOE BLACK	2 PLY NON-HOLE WHITE	2 PLY HOLE-IN-TOE WHITE	2 PLY NON-HOLE BLACK	2 PLY HOLE-IN-TOE BLACK	COMP SIZE	
RK1P-23-10	ORK1P-23-10	BRK1P-23-10	OBRK1P-23-10	RKDF-23-10	ORKDF-23-10	BRKDF-23-10	OBRKDF-23-10	RK ORIGINAL	
RK1P-23-12	ORK1P-23-12	BRK1P-23-12	OBRK1P-23-12	RKDF-23-12	ORKDF-23-12	BRKDF-23-12	OBRKDF-23-12	RK ORIGINAL	
RK1P-23-16	ORK1P-23-16	BRK1P-23-16	OBRK1P-23-16	RKDF-23-16	ORKDF-23-16	BRKDF-23-16	OBRKDF-23-16	RK ORIGINAL	
RK1P-23-20	ORK1P-23-20	BRK1P-23-20	OBRK1P-23-20	RKDF-23-20	ORKDF-23-20	BRKDF-23-20	OBRKDF-23-20	RK ORIGINAL	

35 (X-LARGE)		9"-11" TOP X 5.5"-7" TOE- FLAT MEASUREMENT							
1 PLY NON-HOLE WHITE	1 PLY HOLE-IN-TOE WHITE	1 PLY NON-HOLE BLACK	1 PLY HOLE-IN-TOE BLACK	2 PLY NON-HOLE WHITE	2 PLY HOLE-IN-TOE WHITE	2 PLY NON-HOLE BLACK	2 PLY HOLE-IN-TOE BLACK	COMP SIZE	
RK1P-35-12	ORK1P-35-12	BRK1P-35-12	OBRK1P-35-12	RKDF-35-12	ORKDF-35-12	BRKDF-35-12	OBRKDF-35-12	WIDE SHORT	
RK1P-35-16	ORK1P-35-16	BRK1P-35-16	OBRK1P-35-16	RKDF-35-16	ORKDF-35-16	BRKDF-35-16	OBRKDF-35-16	WIDE MEDIUM	
RK1P-35-20	ORK1P-35-20	BRK1P-35-20	OBRK1P-35-20	RKDF-35-20	ORKDF-35-20	BRKDF-35-20	OBRKDF-35-20	WIDE LONG	

The Royal Flush

Our best hand for softness, cushioning, comfort and fit. This performance fabric is designed to quickly move moisture away from the skin and the breathability of the fiber keeps you drier and more comfortable. With proper care this sock stays soft and resists shrinking, odors, and mildew.

L-CODE: L8430 - AK L-CODE: L8420 - BK

See below for our easy add-on **OPTIONS**.

8 PLY SOCKS ARE NON-RETURNABLE

HEAVY—(Add “H” in front of product code) Ex: (HRF3), (BHRF3), (OBHRF3)

SILVER—(Add “X” in front of product code) Ex: (XRF3), (OXRF3)

- Non-Roll Top
- Form Fitting
- Stretches in length and width
- Round Toe
- Wash and Dry on Low Heat
- Most Sizes fit Symes and PTB Patients



A/B (X-SMALL)		4"-5" TOP X 3"-3.5" TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-A/B-06	ORF*-A/B-06	GORF*-A/B-06	BRF*-A/B-06	OBRF*-A/B-06	GOBRF*-A/B-06	CHILD X-SHORT
RF*-A/B-10	ORF*-A/B-10	GORF*-A/B-10	BRF*-A/B-10	OBRF*-A/B-10	GOBRF*-A/B-10	CHILD SHORT
RF*-A/B-12	ORF*-A/B-12	GORF*-A/B-12	BRF*-A/B-12	OBRF*-A/B-12	GOBRF*-A/B-12	CHILD MEDIUM
RF*-A/B-14	ORF*-A/B-14	GORF*-A/B-14	BRF*-A/B-14	OBRF*-A/B-14	GOBRF*-A/B-14	RK ORIGINAL
RF*-A/B-18	ORF*-A/B-18	GORF*-A/B-18	BRF*-A/B-18	OBRF*-A/B-18	GOBRF*-A/B-18	RK ORIGINAL
RF*-A/B-22	ORF*-A/B-22	GORF*-A/B-22	BRF*-A/B-22	OBRF*-A/B-22	GOBRF*-A/B-22	RK ORIGINAL

B/O (X-SMALL / SMALL)		5"-6" TOP X 3"-3.5" TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-B/0-10	ORF*-B/0-10	GORF*-B/0-10	RF*-B/0-10	OBRF*-B/0-10	GOBRF*-B/0-10	RK ORIGINAL
RF*-B/0-12	ORF*-B/0-12	GORF*-B/0-12	RF*-B/0-12	OBRF*-B/0-12	GOBRF*-B/0-12	RK ORIGINAL
RF*-B/0-14	ORF*-B/0-14	GORF*-B/0-14	RF*-B/0-14	OBRF*-B/0-14	GOBRF*-B/0-14	RK ORIGINAL
RF*-B/0-18	ORF*-B/0-18	GORF*-B/0-18	RF*-B/0-18	OBRF*-B/0-18	GOBRF*-B/0-18	RK ORIGINAL
RF*-B/0-22	ORF*-B/0-22	GORF*-B/0-22	RF*-B/0-22	OBRF*-B/0-22	GOBRF*-B/0-22	RK ORIGINAL

O/1 (SMALL)		6"-7" TOP X 3"-4" TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-0/1-10	ORF*-0/1-10	GORF*-0/1-10	BRF*-0/1-10	OBRF*-0/1-10	GOBRF*-0/1-10	NARROW X-SHORT
RF*-0/1-12	ORF*-0/1-12	GORF*-0/1-12	BRF*-0/1-12	OBRF*-0/1-12	GOBRF*-0/1-12	NARROW SHORT
RF*-0/1-14	ORF*-0/1-14	GORF*-0/1-14	BRF*-0/1-14	OBRF*-0/1-14	GOBRF*-0/1-14	NARROW SHORT
RF*-0/1-18	ORF*-0/1-18	GORF*-0/1-18	BRF*-0/1-18	OBRF*-0/1-18	GOBRF*-0/1-18	NARROW MEDIUM
RF*-0/1-22	ORF*-0/1-22	GORF*-0/1-22	BRF*-0/1-22	OBRF*-0/1-22	GOBRF*-0/1-22	NARROW MEDIUM
RF*-0/1-26	ORF*-0/1-26	GORF*-0/1-26	BRF*-0/1-26	OBRF*-0/1-26	GOBRF*-0/1-26	NARROW LONG

O/2 (SMALL/MEDIUM)		6''-7'' TOP X 4''-5'' TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-0/2-10	ORF*-0/2-10	GORF*-0/2-10	BRF*-0/2-10	OBRF*-0/2-10	GOBRF*-0/2-10	RK ORIGINAL
RF*-0/2-12	ORF*-0/2-12	GORF*-0/2-12	BRF*-0/2-12	OBRF*-0/2-12	GOBRF*-0/2-12	RK ORIGINAL
RF*-0/2-14	ORF*-0/2-14	GORF*-0/2-14	BRF*-0/2-14	OBRF*-0/2-14	GOBRF*-0/2-14	RK ORIGINAL
RF*-0/2-18	ORF*-0/2-18	GORF*-0/2-18	BRF*-0/2-18	OBRF*-0/2-18	GOBRF*-0/2-18	RK ORIGINAL
RF*-0/2-22	ORF*-0/2-22	GORF*-0/2-22	BRF*-0/2-22	OBRF*-0/2-22	GOBRF*-0/2-22	RK ORIGINAL
RF*-0/2-26	ORF*-0/2-26	GORF*-0/2-26	BRF*-0/2-26	OBRF*-0/2-26	GOBRF*-0/2-26	RK ORIGINAL

1/2 (MEDIUM)		7''-8'' TOP X 4''-5'' TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-1/2-10	ORF*-1/2-10	GORF*-1/2-10	BRF*-1/2-10	OBRF*-1/2-10	GOBRF*-1/2-10	REGULAR X-SHORT
RF*-1/2-12	ORF*-1/2-12	GORF*-1/2-12	BRF*-1/2-12	OBRF*-1/2-12	GOBRF*-1/2-12	REGULAR SHORT
RF*-1/2-14	ORF*-1/2-14	GORF*-1/2-14	BRF*-1/2-14	OBRF*-1/2-14	GOBRF*-1/2-14	REGULAR SHORT
RF*-1/2-18	ORF*-1/2-18	GORF*-1/2-18	BRF*-1/2-18	OBRF*-1/2-18	GOBRF*-1/2-18	REGULAR MEDIUM
RF*-1/2-22	ORF*-1/2-22	GORF*-1/2-22	BRF*-1/2-22	OBRF*-1/2-22	GOBRF*-1/2-22	REGULAR MEDIUM
RF*-1/2-26	ORF*-1/2-26	GORF*-1/2-26	BRF*-1/2-26	OBRF*-1/2-26	GOBRF*-1/2-26	REGULAR LONG

2/3 (LARGE)		8''-9'' TOP X 5''-6'' TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-2/3-10	ORF*-2/3-10	GORF*-2/3-10	BRF*-2/3-10	OBRF*-2/3-10	GOBRF*-2/3-10	RK ORIGINAL
RF*-2/3-12	ORF*-2/3-12	GORF*-2/3-12	BRF*-2/3-12	OBRF*-2/3-12	GOBRF*-2/3-12	RK ORIGINAL
RF*-2/3-14	ORF*-2/3-14	GORF*-2/3-14	BRF*-2/3-14	OBRF*-2/3-14	GOBRF*-2/3-14	RK ORIGINAL
RF*-2/3-18	ORF*-2/3-18	GORF*-2/3-18	BRF*-2/3-18	OBRF*-2/3-18	GOBRF*-2/3-18	RK ORIGINAL
RF*-2/3-22	ORF*-2/3-22	GORF*-2/3-22	BRF*-2/3-22	OBRF*-2/3-22	GOBRF*-2/3-22	RK ORIGINAL
RF*-2/3-26	ORF*-2/3-26	GORF*-2/3-26	BRF*-2/3-26	OBRF*-2/3-26	GOBRF*-2/3-26	RK ORIGINAL

3/4 (X-LARGE)		9''-10'' TOP X 6''-7'' TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-3/4-10	ORF*-3/4-10	GORF*-3/4-10	BRF*-3/4-10	OBRF*-3/4-10	GOBRF*-3/4-10	WIDE X-SHORT
RF*-3/4-12	ORF*-3/4-12	GORF*-3/4-12	BRF*-3/4-12	OBRF*-3/4-12	GOBRF*-3/4-12	WIDE SHORT
RF*-3/4-14	ORF*-3/4-14	GORF*-3/4-14	BRF*-3/4-14	OBRF*-3/4-14	GOBRF*-3/4-14	WIDE SHORT
RF*-3/4-18	ORF*-3/4-18	GORF*-3/4-18	BRF*-3/4-18	OBRF*-3/4-18	GOBRF*-3/4-18	WIDE MEDIUM
RF*-3/4-22	ORF*-3/4-22	GORF*-3/4-22	BRF*-3/4-22	OBRF*-3/4-22	GOBRF*-3/4-22	WIDE MEDIUM
RF*-3/4-26	ORF*-3/4-26	GORF*-3/4-26	BRF*-3/4-26	OBRF*-3/4-26	GOBRF*-3/4-26	RK ORIGINAL

4/6 (XX-LARGE)		10''-12'' TOP X 7''-8'' TOE- FLAT MEASUREMENT				* = (3,5,6,8 ply)
NON-HOLE WHITE	HOLE-IN-TOE WHITE	REINFORCED-TOE WHITE	NON-HOLE BLACK	HOLE-IN-TOE BLACK	REINFORCED-TOE BLACK	COMP SIZE
RF*-4/6-10	ORF*-4/6-10	GORF*-4/6-10	BRF*-4/6-10	OBRF*-4/6-10	GOBRF*-4/6-10	X-WIDE X-SHORT
RF*-4/6-12	ORF*-4/6-12	GORF*-4/6-12	BRF*-4/6-12	OBRF*-4/6-12	GOBRF*-4/6-12	X-WIDE SHORT
RF*-4/6-14	ORF*-4/6-14	GORF*-4/6-14	BRF*-4/6-14	OBRF*-4/6-14	GOBRF*-4/6-14	X-WIDE SHORT
RF*-4/6-18	ORF*-4/6-18	GORF*-4/6-18	BRF*-4/6-18	OBRF*-4/6-18	GOBRF*-4/6-18	X-WIDE MEDIUM
RF*-4/6-22	ORF*-4/6-22	GORF*-4/6-22	BRF*-4/6-22	OBRF*-4/6-22	GOBRF*-4/6-22	X-WIDE MEDIUM
RF*-4/6-26	ORF*-4/6-26	GORF*-4/6-26	BRF*-4/6-26	OBRF*-4/6-26	GOBRF*-4/6-26	RK ORIGINAL

THE ROYAL 12-PACK

The Royal Flush product line that you have come to know, and love is now available in a 12-pack! You can now receive all of your patients socks in a single pack. This should be more convenient for you and your patients, and make inventory management easier.

Features of The Royal Flush Sock:

- Comes with 6 (1-ply), 3 (3-ply), 3 (5-ply socks).
- **Non-Roll Top!!** - Eliminates bunching for a smoother, more comfortable fit.
- Available with or without a Distal Opening and/or a Reinforced Distal Opening.
- Available in White or Black
- Form fitting
- Wicks moisture away from the skin.
- Easy care and shrinkage resistant.



6 (1-ply), 3 (3-ply), 3 (5-ply)

Without Hole	With Hole	Size	Fit Sock Size	Length
RF135-01-10	ORF135-01-10	Narrow X-Short	#0, #1	8"-10"
RF135-01-12	ORF135-01-12	Narrow Short	#0, #1	10"-12"
RF135-01-14	ORF135-01-14	Narrow Short+	#0,#1	12"-14"
RF135-12-10	ORF135-12-10	Medium X-Short	#1, #2	8" - 10"
RF135-12-12	ORF135-12-12	Medium Short	#1, #2	10" - 12"
RF135-12-14	ORF135-12-14	Medium Short+	#1, #2	12" - 14"
RF135-23-10	ORF135-23-10	Large X-Short	#2, #3	8" - 10"
RF135-23-12	ORF135-23-12	Large Short	#2, #3	10" - 12"
RF135-23-14	ORF135-23-14	Large Short+	#2, #3	12" - 14"
RF135-34-10	ORF135-34-10	X-Large X-Short	#3 +	8" - 10"
RF135-34-12	ORF135-34-12	X-Large Short	#3 +	10" - 12"
RF135-34-14	ORF135-34-14	X-Large Short+	#3 +	12" - 14"

SIZING CHARTS

6 (2-ply), 3 (3-ply), 3 (5-ply)

Without Hole	With Hole	Size	Fit Sock Size	Length
RF235-01-10	ORF235-01-10	Narrow X-Short	#0, #1	8"-10"
RF235-01-12	ORF235-01-12	Narrow Short	#0, #1	10"-12"
RF235-01-14	ORF235-01-14	Narrow Short+	#0,#1	12"-14"
RF235-12-10	ORF235-12-10	Medium X-Short	#1, #2	8" - 10"
RF235-12-12	ORF235-12-12	Medium Short	#1, #2	10" - 12"
RF235-12-14	ORF235-12-14	Medium Short+	#1, #2	12" - 14"
RF235-23-10	ORF235-23-10	Large X-Short	#2, #3	8" - 10"
RF235-23-12	ORF235-23-12	Large Short	#2, #3	10" - 12"
RF235-23-14	ORF235-23-14	Large Short+	#2, #3	12" - 14"
RF235-34-10	ORF235-34-10	X-Large X-Short	#3 +	8" - 10"
RF235-34-12	ORF235-34-12	X-Large Short	#3 +	10" - 12"
RF235-34-14	ORF235-34-14	X-Large Short+	#3 +	12" - 14"

Add a "H" in front of the "R" to get a Heavyweight
Add a "B" in front of the "R" or the "H" to get a black sock

Royal Open Toe Sock

The Royal Open Toe sock is designed for any application that requires no distal end. May be used to improve the fit of seal-in liners if there is any volume reduction in the proximal area of the residual limb.

Features of the Royal Open Toe Sock

- Non-Roll Top
- Non-Roll Toe Opening.
- Available in 1-6 ply and in 3-6 ply heavy.
- Available in white or black.
- Additional sizes available up request .



White

SIZING CHARTS

1-PLY	2-PLY	3-PLY	5-PLY	6-PLY	Size	Length
SRK1P-01-10	SRKDF-01-10	SRF3-0/1-10	SRF5-0/1-10	SRF6-0/1-10	Narrow X-Short	8"-10"
SRK1P-01-12	SRKDF-01-12	SRF3-0/1-12	SRF5-0/1-12	SRF6-0/1-12	Narrow Short	10"-12"
NA	NA	SRF3-0/1-14	SRF3-0/1-14	SRF6-0/1-14	Narrow Short+	12"-14"
SRK1P-12-10	SRKDF-12-10	SRF3-1/2-10	SRF5-1/2-10	SRF6-1/2-10	Medium X-Short	8" - 10"
SRK1P-12-12	SRKDF-12-12	SRF3-1/2-12	SRF5-1/2-12	SRF6-1/2-12	Medium Short	10" - 12"
NA	NA	SRF3-1/2-14	SRF5-1/2-14	SRF6-1/2-14	Medium Short+	12" - 14"
SRK1P-23-10	SRKDF-23-10	SRF3-2/3-10	SRF5-2/3-10	SRF6-2/3-10	Large X-Short	8" - 10"
SRK1P-23-12	SRKDF-23-12	SRF3-2/3-12	SRF5-2/3-12	SRF6-2/3-12	Large Short	10" - 12"
NA	NA	SRF3-2/3-14	SRF5-2/3-14	SRF6-2/3-14	Large Short+	12" - 14"
SRK1P-35-10	SRKDF-35-10	SRF3-3/4-10	SRF5-3/4-10	SRF6-3/4-10	X-Large X-Short	8" - 10"
SRK1P-35-12	SRKDF-35-12	SRF3-3/4-12	SRF5-3/4-12	SRF6-3/4-12	X-Large Short	10" - 12"
NA	NA	SRF3-3/4-14	SRF5-3/4-14	SRF6-3/4-14	X-Large Short+	12" - 14"

Black

1-PLY	2-PLY	3-PLY	5-PLY	6-PLY	Size	Length
SBRK1P-01-10	SBRKDF-01-10	SBRF3-0/1-10	SBRF5-0/1-10	SBRF6-0/1-10	Narrow X-Short	8"-10"
SBRK1P-01-12	SBRKDF-01-12	SBRF3-0/1-12	SBRF5-0/1-12	SBRF6-0/1-12	Narrow Short	10"-12"
NA	NA	SBRF3-0/1-14	SBRF3-0/1-14	SBRF6-0/1-14	Narrow Short+	12"-14"
SBRK1P-12-10	SBRKDF-12-10	SBRF3-1/2-10	SBRF5-1/2-10	SBRF6-1/2-10	Medium X-Short	8" - 10"
SBRK1P-12-12	SBRKDF-12-12	SBRF3-1/2-12	SBRF5-1/2-12	SBRF6-1/2-12	Medium Short	10" - 12"
NA	NA	SBRF3-1/2-14	SBRF5-1/2-14	SBRF6-1/2-14	Medium Short+	12" - 14"
SBRK1P-23-10	SBRKDF-23-10	SBRF3-2/3-10	SBRF5-2/3-10	SBRF6-2/3-10	Large X-Short	8" - 10"
SBRK1P-23-12	SBRKDF-23-12	SBRF3-2/3-12	SBRF5-2/3-12	SBRF6-2/3-12	Large Short	10" - 12"
NA	NA	SBRF3-2/3-14	SBRF5-2/3-14	SBRF6-2/3-14	Large Short+	12" - 14"
SBRK1P-35-10	SBRKDF-35-10	SBRF3-3/4-10	SBRF5-3/4-10	SBRF6-3/4-10	X-Large X-Short	8" - 10"
SBRK1P-35-12	SBRKDF-35-12	SBRF3-3/4-12	SBRF5-3/4-12	SBRF6-3/4-12	X-Large Short	10" - 12"
NA	NA	SBRF3-3/4-14	SBRF5-3/4-14	SBRF6-3/4-14	X-Large Short+	12" - 14"

THE ROYAL FLUSH MULTI-PLY - 3, 5, 6 PLY COMBINATION SOCK

Multi-Ply is designed to control different volume changes a patient experiences during day to day activities. **Example:** If a patient shrinks distally, instead of adding a short 2 or 3 ply sock, we recommend using a multi-ply sock. The toe is a third of the sock.

Multi-Ply is the same great price and sizing as The Royal Flush sock.

See below for our easy add-on **OPTIONS**

BLACK—(Add “B” in front of product code) Ex: (BRF3/5)

SILVER—(Add “X” in front of product code) Ex: (XRF3/5)


DISTAL OPENING—(Add “O” in front of product code) Ex: (ORF3/5)

REINFORCED DISTAL OPENING—(Add “GO” in front of product code). Ex: (GORF3/5)

EXAMPLE: 5 PLY TOE, 3 PLY TOP, SIZE MEDIUM, 10 TO 12 INCHES LONG.

CODE: RF5/3-1/2-12


12 inch length divided by 3
equals 4 inch toe.



EXAMPLE: 3 PLY TOE, 5 PLY TOP, SIZE MEDIUM, 10 TO 12 INCHES LONG.

CODE: RF3/5-1/2-12

12 inch length divided by 3
equals 4 inch toe.



The Royal Stretch Wool

Our best hand for softness, cushioning, comfort and fit. This performance fabric is designed to quickly move moisture away from the skin and the breathability of the fiber keeps you drier and more comfortable. With proper care this sock stays soft and resists shrinking, odors, and mildew.

L-CODE: L8430 - AK L-CODE: L8420 - BK

See below for our easy add-on **OPTIONS**.

DISTAL OPENING—(Add “O” in front of product code) Ex: (ORS3)

REINFORCED DISTAL OPENING—(Add “GO” in front of product code). Ex: (GORS3)

CUT TOE—(Add “C” in front of product code) Ex: (CRS3)

SEWN TOP—(Add “S” in front of product code) Ex: (SRS3)

- Stretches in Length and Width
- Form Fitting
- Most Sizes fit Symes and PTB Patients
- Round Toe
- Wash and Dry Low Heat



A/B (X-SMALL)		4”-5” TOP X 3”-3.5” TOE- FLAT MEASUREMENT				
3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-A/B-06	ORS3-A/B-06	RS5-A/B-06	ORS5-A/B-06	RS6-A/B-06	ORS6-A/B-06	CHILD X-SHORT
RS3-A/B-10	ORS3-A/B-10	RS5-A/B-10	ORS5-A/B-10	RS6-A/B-10	ORS6-A/B-10	CHILD SHORT
RS3-A/B-12	ORS3-A/B-12	RS5-A/B-12	ORS5-A/B-12	RS6-A/B-12	ORS6-A/B-12	CHILD MEDIUM
RS3-A/B-14	ORS3-A/B-14	RS5-A/B-14	ORS5-A/B-14	RS6-A/B-14	ORS6-A/B-14	RK ORIGINAL
RS3-A/B-18	ORS3-A/B-18	RS5-A/B-18	ORS5-A/B-18	RS6-A/B-18	ORS6-A/B-18	RK ORIGINAL
RS3-A/B-22	ORS3-A/B-22	RS5-A/B-22	ORS5-A/B-22	RS6-A/B-22	ORS6-A/B-22	RK ORIGINAL

B/O (X-SMALL / SMALL)		5”-6” TOP X 3”-3.5” TOE- FLAT MEASUREMENT				
3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-B/0-10	ORS3-B/0-10	RS5-B/0-10	ORS5-B/0-10	RS6-B/0-10	ORS6-B/0-10	RK ORIGINAL
RS3-B/0-12	ORS3-B/0-12	RS5-B/0-12	ORS5-B/0-12	RS6-B/0-12	ORS6-B/0-12	RK ORIGINAL
RS3-B/0-14	ORS3-B/0-14	RS5-B/0-14	ORS5-B/0-14	RS6-B/0-14	ORS6-B/0-14	RK ORIGINAL
RS3-B/0-18	ORS3-B/0-18	RS5-B/0-18	ORS5-B/0-18	RS6-B/0-18	ORS6-B/0-18	RK ORIGINAL
RS3-B/0-22	ORS3-B/0-22	RS5-B/0-22	ORS5-B/0-22	RS6-B/0-22	ORS6-B/0-22	RK ORIGINAL

O/1 (SMALL)		6”-7” TOP X 3”-4” TOE- FLAT MEASUREMENT				
3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-0/1-10	ORS3-0/1-10	RS5-0/1-10	ORS5-0/1-10	RS6-0/1-10	ORS6-0/1-10	NARROW X-SHORT
RS3-0/1-12	ORS3-0/1-12	RS5-0/1-12	ORS5-0/1-12	RS6-0/1-12	ORS6-0/1-12	NARROW SHORT
RS3-0/1-14	ORS3-0/1-14	RS5-0/1-14	ORS5-0/1-14	RS6-0/1-14	ORS6-0/1-14	NARROW SHORT
RS3-0/1-18	ORS3-0/1-18	RS5-0/1-18	ORS5-0/1-18	RS6-0/1-18	ORS6-0/1-18	NARROW MEDIUM
RS3-0/1-22	ORS3-0/1-22	RS5-0/1-22	ORS5-0/1-22	RS6-0/1-22	ORS6-0/1-22	NARROW MEDIUM
RS3-0/1-26	ORS3-0/1-26	RS5-0/1-26	ORS5-0/1-26	RS6-0/1-26	ORS6-0/1-26	NARROW LONG

0/2 (SMALL/MEDIUM)**6"-7" TOP X 4"-5" TOE- FLAT MEASUREMENT**

3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-0/2-10	ORS3-0/2-10	RS5-0/2-10	ORS5-0/2-10	RS6-0/2-10	ORS6-0/2-10	RK ORIGINAL
RS3-0/2-12	ORS3-0/2-12	RS5-0/2-12	ORS5-0/2-12	RS6-0/2-12	ORS6-0/2-12	RK ORIGINAL
RS3-0/2-14	ORS3-0/2-14	RS5-0/2-14	ORS5-0/2-14	RS6-0/2-14	ORS6-0/2-14	RK ORIGINAL
RS3-0/2-18	ORS3-0/2-18	RS5-0/2-18	ORS5-0/2-18	RS6-0/2-18	ORS6-0/2-18	RK ORIGINAL
RS3-0/2-22	ORS3-0/2-22	RS5-0/2-22	ORS5-0/2-22	RS6-0/2-22	ORS6-0/2-22	RK ORIGINAL
RS3-0/2-26	ORS3-0/2-26	RS5-0/2-26	ORS5-0/2-26	RS6-0/2-26	ORS6-0/2-26	RK ORIGINAL

1/2 (MEDIUM)**7"-8" TOP X 4"-5" TOE- FLAT MEASUREMENT**

3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-1/2-10	ORS3-1/2-10	RS5-1/2-10	ORS5-1/2-10	RS6-1/2-10	ORS6-1/2-10	REGULAR X-SHORT
RS3-1/2-12	ORS3-1/2-12	RS5-1/2-12	ORS5-1/2-12	RS6-1/2-12	ORS6-1/2-12	REGULAR SHORT
RS3-1/2-14	ORS3-1/2-14	RS5-1/2-14	ORS5-1/2-14	RS6-1/2-14	ORS6-1/2-14	REGULAR SHORT
RS3-1/2-18	ORS3-1/2-18	RS5-1/2-18	ORS5-1/2-18	RS6-1/2-18	ORS6-1/2-18	REGULAR MEDIUM
RS3-1/2-22	ORS3-1/2-22	RS5-1/2-22	ORS5-1/2-22	RS6-1/2-22	ORS6-1/2-22	REGULAR MEDIUM
RS3-1/2-26	ORS3-1/2-26	RS5-1/2-26	ORS5-1/2-26	RS6-1/2-26	ORS6-1/2-26	REGULAR LONG

2/3 (LARGE)**8"-9" TOP X 5"-6" TOE- FLAT MEASUREMENT**

3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-2/3-10	ORS3-2/3-10	RS5-2/3-10	ORS5-2/3-10	RS6-2/3-10	ORS6-2/3-10	RK ORIGINAL
RS3-2/3-12	ORS3-2/3-12	RS5-2/3-12	ORS5-2/3-12	RS6-2/3-12	ORS6-2/3-12	RK ORIGINAL
RS3-2/3-14	ORS3-2/3-14	RS5-2/3-14	ORS5-2/3-14	RS6-2/3-14	ORS6-2/3-14	RK ORIGINAL
RS3-2/3-18	ORS3-2/3-18	RS5-2/3-18	ORS5-2/3-18	RS6-2/3-18	ORS6-2/3-18	RK ORIGINAL
RS3-2/3-22	ORS3-2/3-22	RS5-2/3-22	ORS5-2/3-22	RS6-2/3-22	ORS6-2/3-22	RK ORIGINAL
RS3-2/3-26	ORS3-2/3-26	RS5-2/3-26	ORS5-2/3-26	RS6-2/3-26	ORS6-2/3-26	RK ORIGINAL

3/4 (X-LARGE)**9"-10" TOP X 6"-7" TOE- FLAT MEASUREMENT**

3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-3/4-10	ORS3-3/4-10	RS5-3/4-10	ORS5-3/4-10	RS6-3/4-10	ORS6-3/4-10	WIDE X-SHORT
RS3-3/4-12	ORS3-3/4-12	RS5-3/4-12	ORS5-3/4-12	RS6-3/4-12	ORS6-3/4-12	WIDE SHORT
RS3-3/4-14	ORS3-3/4-14	RS5-3/4-14	ORS5-3/4-14	RS6-3/4-14	ORS6-3/4-14	WIDE SHORT
RS3-3/4-18	ORS3-3/4-18	RS5-3/4-18	ORS5-3/4-18	RS6-3/4-18	ORS6-3/4-18	WIDE MEDIUM
RS3-3/4-22	ORS3-3/4-22	RS5-3/4-22	ORS5-3/4-22	RS6-3/4-22	ORS6-3/4-22	WIDE MEDIUM
RS3-3/4-26	ORS3-3/4-26	RS5-3/4-26	ORS5-3/4-26	RS6-3/4-26	ORS6-3/4-26	RK ORIGINAL

4/6 (XX-LARGE)**10"-12" TOP X 7"-8" TOE- FLAT MEASUREMENT**

3 PLY NON-HOLE	3 PLY HOLE-IN-TOE	5 PLY NON-HOLE	5 PLY HOLE-IN-TOE	6 PLY NON-HOLE	6 PLY HOLE-IN-TOE	COMP SIZE
RS3-4/6-10	ORS3-4/6-10	RS5-4/6-10	ORS5-4/6-10	RS6-4/6-10	ORS6-4/6-10	X-WIDE X-SHORT
RS3-4/6-12	ORS3-4/6-12	RS5-4/6-12	ORS5-4/6-12	RS6-4/6-12	ORS6-4/6-12	X-WIDE SHORT
RS3-4/6-14	ORS3-4/6-14	RS5-4/6-14	ORS5-4/6-14	RS6-4/6-14	ORS6-4/6-14	X-WIDE SHORT
RS3-4/6-18	ORS3-4/6-18	RS5-4/6-18	ORS5-4/6-18	RS6-4/6-18	ORS6-4/6-18	X-WIDE MEDIUM
RS3-4/6-22	ORS3-4/6-22	RS5-4/6-22	ORS5-4/6-22	RS6-4/6-22	ORS6-4/6-22	X-WIDE MEDIUM
RS3-4/6-26	ORS3-4/6-26	RS5-4/6-26	ORS5-4/6-26	RS6-4/6-26	ORS6-4/6-26	RK ORIGINAL

POST OP PROSTHETIC SOCK

Full Fashioned Post Op Prosthetic Sock

Designed for immediate post operative surgery. This sock has a seamless round toe, versus the tradition square toe, engineered for a more conforming fit for the post operative patient. (Please indicated ROLLED if this is desired)

See below for our easy add-on **OPTIONS**

BLACK—(Add “B” in front of product code) Ex: **(BRF3)**

SILVER—(Add “X” in front of product code) Ex: **(XRF3)**



ALL MEASUREMENTS LYING FLAT

CODE	TOP x TOE	DESCRIPTION
RF3-00-14	5" x 3.5"	(12"-14") Royal Post Op 3 Ply White X-Small No Hole
RF3-00-18	5" x 4"	(16"-18") Royal Post Op 3 Ply White Small No Hole
RF3-10-18	6" x 5"	(16"-18") Royal Post Op 3 Ply White Medium No Hole
RF3-20-18	8" x 6"	(16"-18") Royal Post Op 3 Ply White Large No Hole
RF3-30-18	9" x 7"	(16"-18") Royal Post Op 3 Ply White X-Large No Hole

ALL MEASUREMENTS LYING FLAT

CODE	TOP x TOE	DESCRIPTION
RF5-00-14	5" x 3.5"	(12"-14") Royal Post Op 5 Ply White X-Small No Hole
RF5-00-18	5" x 4"	(16"-18") Royal Post Op 5 Ply White Small No Hole
RF5-10-18	6" x 5"	(16"-18") Royal Post Op 5 Ply White Medium No Hole
RF5-20-18	8" x 6"	(16"-18") Royal Post Op 5 Ply White Large No Hole
RF5-30-18	9" x 7"	(16"-18") Royal Post Op 5 Ply White X-Large No Hole



HEMIPELVECTOMY/HIP DISARTICULATION SOCK

The Royal Knit Full Fashioned Hip Sock

Utilizing a revolutionary knitting technique, we have developed a hip disarticulation sock which makes fitting and ordering easy as possible.

This product features fully seamless knitting in four sizes, will fit ninety-five percent of your patients. This includes right and left amputation, both male and female.

The only measurement required to order this product is waist circumference. If a custodial opening is desired (add \$4.00 to price below), simply designate right or left amputation. Bi-lateral also known as a bucket sock available on request.

Available exclusively in 3-ply. This product features excellent wicking action, quick drying, overall softness and a more conforming fit.

See below for our easy add-on **OPTIONS**

BLACK-(Add “B” in front of the product code) Ex: **(BRF3-HP)**

CUSTODIAL OPENING-(Add “F” in front of product code) Ex: **(FRF3-HP)**

BI-LATERAL-(Add “BL” in the place of the “HP” in the product code) Ex: **(RF3-BL)**

Please refer to the sizing chart below to place your order.

HIP SOCK			
CODE	WAIST SIZE CIRC.	LENGTH	DESCRIPTION
RF3-HP-XS	15”-21”	12”-14”	Royal Hip Sock 3 Ply White X-Small No Fly
RF3-HP-SM	22”-28”	14”-16”	Royal Hip Sock 3 Ply White Small No Fly
RF3-HP-MD	28”-36”	16”-18”	Royal Hip Sock 3 Ply White Medium No Fly
RF3-HP-LG	36”-47”	16”-18”	Royal Hip Sock 3 Ply White Large No Fly



Torso Body Sock

"Fashion By Design"

The full fashion body sock can be ordered in our standard design with or without shoulder caps. The body sock can also be ordered with options that include V-neck, Tube Top and axilla flaps.

Full fashioned seamless knitting of the Royal Knit body sock means no unraveling stitches. The excellent wicking action of the material quickly moves moisture away from the skin to the outer layer of the fabric. The breathability of the fiber provides a drier, softer and more comfortable fit. Consequent to the wicking action created by Coolmax® the possibility of contact dermatitis is greatly reduced.

Only five overlapping sizes are available for convenient ordering.

Machine wash and tumble dry warm. The Royal Knit body sock will accept normal bleach applications. Sheet fabric softener recommended. **DO NOT HAND WASH.**

To order, please refer to the body sock sizing chart on next page.



V-Neck
(WITH OR WITHOUT
AXILLA FLAPS)



Crew Neck
(WITH OR WITHOUT
AXILLA FLAPS)

SIZING CHART

L-CODE: L0984

CREW NECK

WHITE	SILVER	WIDTH	LENGTH	FITS CIRCUMFERENCE
TLSO-XS-12	XTLSO-XS-12	4" (10cm)	12" (30cm)	12"-20" (30-51cm)
TLSO-SM-14	XTLSO-SM-14	6" (15cm)	14" (35cm)	16"-27" (41-68cm)
TLSO-MD-16	XTLSO-MD-16	8" (20cm)	16" (41cm)	22"-37" (56-94cm)
TLSO-ML-20	XTLSO-ML-20	8" (20cm)	20" (51cm)	22"-37" (56-94cm)
TLSO-LG-20	XTLSO-LG-20	10" (25cm)	20" (51cm)	30"-50" (76-127cm)
TLSO-LL-24	XTLSO-LL-24	10" (25cm)	24" (61cm)	30"-50" (76-127cm)
TLSO-XL-24	XTLSO-SL-24	12" (30cm)	24" (61cm)	36"-60" (91-165cm)

CREW NECK (AXILLA FLAPS)

WHITE	SILVER	WIDTH	LENGTH	FITS CIRCUMFERENCE
FTLSO-XS-12	FXTLSO-XS-12	4" (10cm)	12" (30cm)	12"-20" (30-51cm)
FTLSO-SM-14	FXTLSO-SM-14	6" (15cm)	14" (35cm)	16"-27" (41-68cm)
FTLSO-MD-16	FXTLSO-MD-16	8" (20cm)	16" (41cm)	22"-37" (56-94cm)
FTLSO-ML-20	FXTLSO-ML-20	8" (20cm)	20" (51cm)	22"-37" (56-94cm)
FTLSO-LG-20	FXTLSO-LG-20	10" (25cm)	20" (51cm)	30"-50" (76-127cm)
FTLSO-LL-24	FXTLSO-LL-24	10" (25cm)	24" (61cm)	30"-50" (76-127cm)
FTLSO-XL-24	FXTLSO-SL-24	12" (30cm)	24" (61cm)	36"-60" (91-165cm)


V-NECK

WHITE	SILVER	WIDTH	LENGTH	FITS CIRCUMFERENCE
VTLSO-XS-12	VXTLSO-XS-12	4" (10cm)	12" (30cm)	12"-20" (30-51cm)
VTLSO-SM-14	VXTLSO-SM-14	6" (15cm)	14" (35cm)	16"-27" (41-68cm)
VTLSO-MD-16	VXTLSO-MD-16	8" (20cm)	16" (41cm)	22"-37" (56-94cm)
VTLSO-ML-20	VXTLSO-ML-20	8" (20cm)	20" (51cm)	22"-37" (56-94cm)
VTLSO-LG-20	VXTLSO-LG-20	10" (25cm)	20" (51cm)	30"-50" (76-127cm)
VTLSO-LL-24	VXTLSO-LL-24	10" (25cm)	24" (61cm)	30"-50" (76-127cm)
VTLSO-XL-24	VXTLSO-SL-24	12" (30cm)	24" (61cm)	36"-60" (91-165cm)

V-NECK (AXILLA FLAPS)

WHITE	SILVER	WIDTH	LENGTH	FITS CIRCUMFERENCE
VFTLSO-XS-12	VFXTLSO-XS-12	4" (10cm)	12" (30cm)	12"-20" (30-51cm)
VFTLSO-SM-14	VFXTLSO-SM-14	6" (15cm)	14" (35cm)	16"-27" (41-68cm)
VFTLSO-MD-16	VFXTLSO-MD-16	8" (20cm)	16" (41cm)	22"-37" (56-94cm)
VFTLSO-ML-20	VFXTLSO-ML-20	8" (20cm)	20" (51cm)	22"-37" (56-94cm)
VFTLSO-LG-20	VFXTLSO-LG-20	10" (25cm)	20" (51cm)	30"-50" (76-127cm)
VFTLSO-LL-24	VFXTLSO-LL-24	10" (25cm)	24" (61cm)	30"-50" (76-127cm)
VFTLSO-XL-24	VFXTLSO-SL-24	12" (30cm)	24" (61cm)	36"-60" (91-165cm)

The Royal Orthotic Brace/Cast Sock



AFO/KAFO

Produced with our full fashioned, seamless knitting techniques in two ply, our brace sock has the same inherent qualities of our standard prosthetic sock products.

Our brace sock allows for a more snug, wrinkle free fit in both AFO and KAFO sizing. Available in open and closed toe design. Machine wash and tumble dry warm. Sheet fabric softener recommended.

See below for our easy add-on **OPTIONS**

SILVER-(Add "X" in front of product code) Ex: **(XROBS)**

OPEN TOE-(Add "S" in front of product code) Ex: **(SROBS)**

ALL MEASUREMENTS LYING FLAT

WHITE	BLACK	TOP x TOE	DESCRIPTION
ROBS-0A-14	BROBS-0A-14	4"-5" x 2.5"-3"	(12"-16") Royal Brace Sock 2 Ply X-Small Closed Toe
ROBS-01-16	BROBS-01-16	5.5"-7" x 3"-4"	(16"-20") Royal Brace Sock 2 Ply Small Closed Toe
ROBS-01-20	BROBS-01-20	5.5"-7" x 3"-4"	(22"-28") Royal Brace Sock 2 Ply Small Closed Toe
ROBS-01-24	BROBS-01-24	5.5"-7" x 3"-4"	(28"-36") Royal Brace Sock 2 Ply Small Closed Toe
ROBS-12-16	BROBS-12-16	7"-8" x 4"-5"	(16"-20") Royal Brace Sock 2 Ply Medium Closed Toe
ROBS-12-20	BROBS-12-20	7"-8" x 4"-5"	(22"-28") Royal Brace Sock 2 Ply Medium Closed Toe
ROBS-12-28	BROBS-12-28	7"-8" x 4"-5"	(32"-40") Royal Brace Sock 2 Ply Medium Closed Toe
ROBS-25-32	BROBS-25-32	MEDIUM/LARGE 36"-44" 2 PLY C/M BRACE SOCK	

FULL FASHIONED FRACTURE/BRACE SOCK

Produced with our full fashioned, seamless, knitting techniques in 3-ply or 5-ply, our fracture cast sock has the same outstanding qualities as our brace sock.

The stretch of our cast sock provides a more conforming, wrinkle free fit in all sizes. The product can be ordered with open or closed toe design and can be used with most plastic and plaster casting materials.

Can be gas autoclaved. Machine wash and tumble dry warm. The Royal Knit cast sock will accept normal bleach applications (white only). Sheet fabric softener recommended.

ALL MEASUREMENTS LYING FLAT

CODE	TOP x TOE	DESCRIPTION
RF3-06-36	6" X 3"	(34"-36") Royal Fracture Cast 3 Ply 36X6 Closed Toe
RF3-09-35	9" X 3.5"	(33"-35") Royal Fracture Cast 3 Ply 35X9 Closed Toe
RF3-09-39	9" X 3.5"	(37"-39") Royal Fracture Cast 3 Ply 39X9 Closed Toe
RF3-09-42	9" X 3.5"	(40"-42") Royal Fracture Cast 3 Ply 42X9 Closed Toe
RF3-09-46	9" X 3.5"	(44"-46") Royal Fracture Cast 3 Ply 46X9 Closed Toe
RF3-10-46	10" X 3.5"	(44"-46") Royal Fracture Cast 3 Ply 46X10 Closed Toe

ALL MEASUREMENTS LYING FLAT

CODE	TOP x TOE	DESCRIPTION
RF5-06-36	6" X 3"	(34"-36") Royal Fracture Cast 5 Ply 36X6 Closed Toe
RF5-09-35	9" X 3.5"	(33"-35") Royal Fracture Cast 5 Ply 35X9 Closed Toe
RF5-09-39	9" X 3.5"	(37"-39") Royal Fracture Cast 5 Ply 36X9 Closed Toe
RF5-09-42	9" X 3.5"	(40"-42") Royal Fracture Cast 5 Ply 42X9 Closed Toe
RF5-09-46	9" X 3.5"	(44"-46") Royal Fracture Cast 5 Ply 46X9 Closed Toe
RF5-10-46	10" X 3.5"	(44"-46") Royal Fracture Cast 5 Ply 46X10 Closed Toe

DIABETIC SOCKS

DIABETIC SOCKS

CODE	DESCRIPTION
DIE-07	(5"-7") Royal Diabetic Socks 2 Ply Medium Closed Toe (Pair)
DIE-09	(7"-9") Royal Diabetic Socks 2 Ply Medium Closed Toe (Pair)

CAN ALSO BE USED AS A PARTIAL FOOT SOCK

AVAILABLE IN BLACK, WHITE & SILVER

CRANIAL CAPS

CRANIAL CAP

CODE	DESCRIPTION
WRKDF-23-10	(8"-12") Double Fleece 2 Ply White Cranial Cap
WRK1P-23-10	(8"-12") Double Fleece 1 Ply White Cranial Cap

FRACTURE INTERFACES

ARM SLEEVE, HUMERAL AND ULNAR SLEEVES



For interfacing casts or braces, the arm sleeve features full fashioned, seamless knitting with or without thumb in two ply.

Only two sizes, children's and adult are required to fit all patients. The excellent wicking action quickly moves moisture away from the skin to the outer layer of the fabric. The breathability of the fiber provides a drier, softer and more comfortable fit.

Machine wash and tumble dry warm. Sheet fabric softener recommended. (white only)

HAND SLEEVES WITH KNITTED THUMB

CODE	DESCRIPTION
WTK2-SM-10	(8"-10") Royal Hand Sleeves 2 Ply White Small Knitted Thumb
WTK2-SM-12	(10"-12") Royal Hand Sleeves 2 Ply White Small Knitted Thumb
WTK2-MD-04	(2"-4") Royal Hand Sleeves 2 Ply White Medium Knitted Thumb
WTK2-MD-10	(8"-10") Royal Hand Sleeves 2 Ply White Medium Knitted Thumb
WTK2-MD-18	(16"-18") Royal Hand Sleeves 2 Ply White Medium Knitted Thumb

HAND SLEEVES WITH THUMB HOLE

CODE	DESCRIPTION
WTH2-SM-10	(8"-10") Royal Hand Sleeves 2 Ply White Small Thumb Hole
WTH2-SM-12	(10"-12") Royal Hand Sleeves 2 Ply White Small Thumb Hole
WTH2-MD-04	(2"-4") Royal Hand Sleeves 2 Ply White Medium Thumb Hole
WTH2-MD-10	(8"-10") Royal Hand Sleeves 2 Ply White Medium Thumb Hole
WTH2-MD-18	(16"-18") Royal Hand Sleeves 2 Ply White Medium Thumb Hole

ARM SLEEVES

CODE	DESCRIPTION
ARM-SM-06	(4"-6") Royal Hand Sleeves 2 Ply White Small Sleeve
ARM-SM-08	(6"-8") Royal Hand Sleeves 2 Ply White Small Sleeve
ARM-MD-10	(8"-10") Royal Hand Sleeves 2 Ply White Medium Sleeve
ARM-MD-12	(10"-12") Royal Hand Sleeves 2 Ply White Medium Sleeve

Hip Orthosis Sock

HIP ORTHOSIS SOCK

This garment fits most adults and is useful with every type of hip orthosis. The stretch allows the longer leg to be positioned first and the garment to be extended and relaxed carefully over sensitive hip areas.

The acrylic fiber tends to move moisture away from the skin surface and the 1-ply fabric is enough to protect the skin from brace parts.

It can be easily machine washed, tumbled dry.

The design provides adequate openings so that no changing is required during the day.

ONE SIZE FITS ALL

PRODUCT CODE #HL1P-11-34

L-CODE 2850



THE Royal Toe Sock

TOE SOCK

The Royal Toe sock features a separated toe to accommodate AFO toe straps. It provides maximum protection from wrinkle and seam pressure by using completely seamless knitting.

Sized for AFO applications. The stretch no-heel design accommodates varying foot length to leg situations.

It can be easily machine washed, tumbled dry.

The material is the highest wicking and softest environment possible.

AVAILABLE IN 3 SIZES

**CHILDS SMALL/REGULAR
TOE-SM-12**

**MEDIUM
TOE-MD-16**

**ADULT SMALL/LARGE
TOE-LG-16**



L-CODE INDEX

For customer convenience we have included the L-codes listed below that correlate to most of our products. We feel obligated to remind our facilities that it is their responsibility to determine the applicability of these codes when using them for the purpose of reimbursement.

L0984 – Protective Body Sock

L2840 – Addition to Lower Extremity Orthosis, Tibial Length Sock

L2850 – Addition to Lower Extremity Orthosis, Femoral Length Sock

L3995 – Upper Extremity Fracture Orthosis, Sock Fracture or Equal

L8420 – Prosthetic Sock, Wool, Below Knee

L8430 – Prosthetic Sock, Wool, Above Knee

L8435 – Prosthetic Sock, Wool, Upper Limb

L8400 – Prosthetic Sheath, Below Knee

L8410 – Prosthetic Sheath, Above Knee

NOTES